


QUICK GUIDE IMPRESSIONS SDIS TL

Color coding		
	blue	Endosteal implant diameter 2.9 mm
	yellow	Endosteal implant diameter 3.3 mm
	red	Endosteal implant diameter 4.1 mm
	green	Endosteal implant diameter 4.8 mm

NNC SP IMPLANT			RN SP IMPLANT			WN SP IMPLANT				
● SP NNC Ø 3.3 mm			● SP RN Ø 3.3 mm ● SP RN Ø 4.1 mm ● SP RN Ø 4.8 mm			● SP WN Ø 4.8 mm				
CONNECTION TYPES NNC Ø 3.5 mm			CONNECTION TYPES RN Ø 4.8 mm			CONNECTION TYPES WN Ø 6.5 mm				
										
Digital Impression	Open Tray	Closed Tray	Digital Impression	Open Tray	Closed Tray	Digital Impression	Open Tray	Closed Tray		
Parts needed 048.383	Parts needed 048.128	Parts needed 048.134	Parts needed 048.384	Parts needed 048.090	Parts needed 048.862	Parts needed 048.385	Parts needed 048.091	Parts needed 048.870		
	 048.127 	 048.127 	 048.129 	 048.010  048.124 	 048.070V4 / 048.017V4   048.124 	 048.172 	 048.171 	 048.095 / 048.013   048.171 		
			ABUTMENT LEVEL			ABUTMENT LEVEL				
			Open Tray		Closed Tray		Open Tray		Closed Tray	
			048.866	048.108	048.864	048.108	048.869	048.872	048.871	048.872
										

QUICK GUIDE IMPRESSIONS SDIS BL/BLT





Color coding		
● blue	Endosteal implant diameter 2.9 mm	
● yellow	Endosteal implant diameter 3.3 mm	
● red	Endosteal implant diameter 4.1 mm	
● green	Endosteal implant diameter 4.8 mm	

SC IMPLANT (BLT)			NC IMPLANT			RC IMPLANT		
CONNECTION TYPES			CONNECTION TYPES			CONNECTION TYPES		
● BLT SC Ø 2.9 mm			● BL NC Ø 3.3 mm ● BLT NC Ø 3.3 mm			● BL RC 4.1 mm ● BLT RC 4.1 mm ● BL RC 4.8 mm ● BLT RC 4.8 mm		
SC Ø 2.9 mm			NC Ø 3.3 mm			RC Ø 4.1 mm RC Ø 4.8 mm		
Digital Impression	Open Tray	Closed Tray	Digital Impression	Open Tray	Closed Tray	Digital Impression	Open Tray	Closed Tray
Parts needed	Parts needed	Parts needed	Parts needed	Parts needed	Parts needed	Parts needed	Parts needed	Parts needed
	025.0021	025.0062	025.0078	025.2202 025.2205	025.2201	025.0079	025.4202 025.4205	025.4201
				 			 	
025.0077	025.0022	025.0020		025.0057* 025.0058*			025.0060* 025.0061*	
				 			 	
025.0024	025.0023	025.0023	025.2102	025.2101	025.2101	025.4102	025.4101	025.4101
								
ABUTMENT LEVEL SRA 3.5 MM			ABUTMENT LEVEL SRA 4.6 MM					
Digital Impression	Open Tray	Closed Tray	Digital Impression	Open Tray	Closed Tray			
025.0080	025.2243	025.2245	025.0081	025.2244	025.2246			
	025.0011*	025.0013*		025.0012*	025.0014*			
025.0007	023.2754	023.2754	025.0008	023.4756	023.4756			

*non engaging

QUICK GUIDE

T L C & T L X

NT Implant TLC Ø 3.3 mm TLC/TLX Ø 3.75 mm TLX Ø 4.5 mm			RT Implant TLC Ø 3.3 mm TLC/TLX Ø 3.75 mm TLX Ø 4.5 mm			WT Implant TLC Ø 4.5 mm TLC/TLX Ø 5.5 mm TLX Ø 6.5 mm		
NT (Ø 3.5 mm)			RT (Ø 4.8 mm)			WT (Ø 6.5 mm)		
Digital impression	Open-tray impression	Closed-tray impression	Digital impression	Open-tray impression	Closed-tray impression	Digital impression	Open-tray impression	Closed-tray impression
Scanbody 036.3230 	For crown 036.0000 036.0001 	For crown 036.0002 	Scanbody 036.3230 	For crown 036.1000 036.1001 	For crown 036.1002 	Scanbody 036.3230 	For crown 036.2000 036.2001 	For crown 036.2002 
	For bridge 065.0140 065.0137 	For bridge 065.0143 		For bridge 065.0141 065.0138 	For bridge 065.0144 		For bridge 065.0142 065.0139 	For bridge 065.0145 
Repositionable analog 036.1102 	036.0100 	036.0100 	Repositionable analog 036.2100 	036.0102 	036.0102 	Repositionable analog 036.2102 	036.1100 	036.1100 




International Headquarters
Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com

© Institut Straumann AG, 2024. All rights reserved.
Straumann® and /or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and /or its affiliates.

450.019/en/D/00 03/25

QUICK GUIDE

B L C & B L X

Implant Level Impression Restoration with RB/WB components BLC Ø3.3mm, 3.75mm, 4.5mm, 5.5mm, 6.5mm BLX Ø 3.5mm, 4.0mm, 4.5mm, 5.0mm, 5.5mm, 6.5mm			Implant Level Impression Restoration with WB components (for wide emergence profile) BLC Ø4.5mm, 5.5mm, 6.5mm BLX Ø5.0mm, 5.5mm, 6.5mm			Abutment Level Impression Screw-Retained Abutments		
Digital impression	Open-tray impression	Closed-tray impression	Digital impression	Open-tray impression	Closed-tray impression	Digital impression	Open-tray impression	Closed-tray impression
Scanbody 065.0103 	For crown 065.0031 065.0033 	For crown 065.0146 065.0148 	Scanbody 065.0103 	For crown 065.0032 065.0034 	For crown 065.4810 	Scanbody 025.0081 	For crown 025.2244 	For crown 025.2246 
Repositionable analog 065.0023 	065.0021 	065.0021 	Repositionable analog 065.0024 	065.0022 	065.0022 	Repositionable analog 025.0008 	For Multi-Unit 025.0012 	For Multi-Unit 025.0014 
							023.4756 025.0050 023.4757 	023.4756 025.0050 023.4757 

450.019/en/D/00 03/25